



Provided for your reference by
cpapXchange.com

AirSense™ 10
AirCurve™ 10

ENGLISH

This Air10 device is not equipped with a wireless module. This addendum replaces content found in the User Guide, Clinical Guide and Quick Start Guide for the following devices:

- AirSense 10 AutoSet
- AirCurve 10 VAuto / ASV / ST / ST-A

User guide

The following information replaces content found in the User Guide.

About the control panel

The following icons are not applicable to your device.



Wireless signal strength (green)



Wireless transfer not enabled (grey)



No wireless connection



Airplane mode

Therapy data

Your device records your therapy data for you and your care provider so they can view and make changes to your therapy if required. The data is recorded and then transferred to your care provider via an SD card.

Data transmission

Wireless data transmission is not available for your device. To transfer your therapy data to your care provider, see the SD card section of the User Guide.

Travelling by plane

When using the device on a plane:

- Airplane mode is not available or required for your device.

Troubleshooting

Problem/possible cause	Solution
My therapy data has not been sent to my care provider	
Wireless data transmission is not available for your device.	Talk to your care provider about transferring therapy data. See the SD card section of the User Guide.

Technical specifications

Wireless module None	FCC / IC ID: The device is not equipped with a wireless module and hence the FCC and ISSED rules are not applicable
-------------------------	---

Clinical guide

The following information replaces content found in the Clinical Guide.

About the control panel

The following icons are not applicable to this device.



Wireless signal strength (green)



Wireless transfer not enabled (grey)



No wireless connection



Airplane mode

Data management and therapy compliance

For therapy management, the device stores patient therapy data on the SD card. This data can be transferred via an SD Card Reader to ResMed's ResScan™ therapy management system. For more information on therapy management with ResScan, refer to the manual supplied with the software.

Remote monitoring

Remote monitoring is not available as this device is not equipped with a wireless module.

Data storage

Data can be transferred to therapy management software via SD card. The different ways of transferring data are detailed in the table below.

Type of Data	Data transfer method		Sessions stored
	SD card to ResScan	SD Card to AirView (card-to-cloud)	
Summary data (compliance data)	✓	✓	365
Night profile data	✓	✓	Limited by usage and SD card storage capacity
High resolution flow and pressure data (25 Hz - every 40 ms)	✓		

Therapy data

Wireless data transmission is not available as this device is not equipped with a wireless module. To transfer the patient's therapy data, see the **SD card** section of the User Guide.

Travelling by plane

When using the device on a plane:

- Airplane mode is not available or required for this device.

Troubleshooting

Problem/possible cause	Solution
The patient's therapy data has not been transmitted	
Wireless data transmission is not available for this device.	Talk to your patient about transferring therapy data. See the SD card section of the User Guide.

Technical specifications

Wireless module None	FCC / IC ID: The device is not equipped with a wireless module and hence the FCC and ISED rules are not applicable
-------------------------	--

Quick start guide

The following information replaces content found in the Quick Start Guide.


Register today

myAir is not supported for your device. Product support is available at ResMed.com/support.

 **ResMed Pty Ltd**
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. Air 10, AirCurve, AirSense, AirView, ResScan and myAir are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip. © 2022 ResMed. 3781046/1 2022-04

ResMed.com



3781046



ResMed

Provided for your reference by
cpapXchange.com

AirSense™ 10 **AUTOSET**

AUTOSET FOR HER

ELITE

CPAP



Home button for returning to the
Home screen at any time

Start/Stop button for starting therapy

Dial for navigating and changing options

Air outlet for connecting the air tubing

Water tub for humidification

Quick start guide
English

Rx Only

Read your entire AirSense™ 10 User Guide before use.

Setup



Place the device on a stable level surface.



Plug the power connector into the rear of the device. Connect one end of the power cord into the power supply unit and the other end into the power outlet.




Connect the air tubing firmly to the air outlet located on the rear of the device.



Open the water tub and fill with water up to the maximum water level mark.



Close the water tub and insert it into the side of the device.



Connect the free end of the air tubing firmly onto the assembled mask. Press Start/Stop to begin therapy.

Adjusting to therapy


If you are a first time user, you might need some time to get used to therapy. This is not unusual as it takes most patients between one and two weeks to adjust to the air pressure.

Overcoming symptoms by changing comfort settings


If you are having difficulties getting used to therapy, use the AirSense 10 comfort features to help you on your journey to better sleep.

- **Dry or runny nose**—If you are getting a dry or runny nose, adjust the Humidity Level by turning it up.
- **Droplets of water (condensation)**—If you are getting droplets of water on your nose, mask or air tubing, adjust the Humidity Level by turning it down.

To change the Humidity Level:




1. Press the dial to enter My Options.



2. Turn the dial to highlight Humidity Level and then press to select it.




3. Turn the dial to adjust to your preferred setting.




4. Press the dial to save the change.

- **Trouble falling asleep**—If you are having trouble falling asleep with high pressure, turn on Auto Ramp or increase Ramp Time.
- **Bloated feeling**—If you are experiencing a slightly bloated feeling from swallowing air, turn on Auto Ramp or increase Ramp Time.
- **Feeling of not getting enough air**—If you feel like you are not getting enough air, turn Ramp Time to Off.

To change the Ramp Time:




1. Press the dial to enter My Options.



2. Turn the dial to highlight Ramp Time and then press to select it.



3. Turn the dial to adjust to your preferred setting.



4. Press the dial to save the change.



Cleaning

1. Wash the water tub and air tubing in warm water using mild detergent.
Do not wash in a dishwasher or washing machine.
2. Rinse the water tub and air tubing thoroughly and allow to dry out of direct sunlight and/or heat.
3. Wipe the exterior of the device with a dry cloth.

Refer to your mask User Guide for detailed instructions on cleaning your mask.

Checking and replacing your parts regularly

It is important for your comfort and health that you check and replace your parts and supplies regularly. Replacing your parts and supplies on a regular basis helps ensure you are receiving optimal therapy and continued comfort.

Refer to your AirSense 10 User Guide for detailed instructions on checking your device. Contact your care provider for information on your replacement schedule.

Replacement schedule guide:

Every Month	Every 3 Months	Every 6 Months
Mask cushion or pillow	Mask frame (excludes headgear)	Headgear
Air filters	Air tubing	Water tub

Register today

myAir™ allows you to access your therapy data and provides personalized support.

Every journey starts with a first step. Register today at ResMed.com/myAir

Serial number	Device number
Mask type and size	Date of setup
Notes	

 **ResMed Ltd**
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

DISTRIBUTED BY
ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA
 ResMed (UK) Ltd 96 Jubilee Ave Milton Park Abingdon Oxfordshire OX14 4RW UK

See www.resmed.com for other ResMed locations worldwide. For patent information, see www.resmed.com/ip.
AirSense and myAir are trademarks of ResMed Ltd. © 2014 ResMed Ltd. 378310/3 2014-07

ResMed.com



378310