

**CPAPXCHANGE**

SLEEP WELL TONIGHT



**React Health's best selling full face mask offers the following features for patient comfort and compliance:**

- No forehead brace allows for increased line of sight - even with your glasses on
- Self-adjusting ergonomic cushion design allows for superior seal
- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- 360-degree dual-swivel elbow for active sleepers and allows for sleep in any position
- Air vent diffuser for a quiet night sleep
- Anti-phyxia valve safety feature opens for fresh air in case of machine shut off



REACTHEALTH

## **Siesta** Full Face Mask



# The Siesta Full Face Mask is labeled on the back of the mask.

Allows for ease of identification for troubleshooting and resupply



## Mask Part Numbers

- SFF1000 - Siesta Full Face Mask - All Sizes Kit
- SFF1001 - Siesta Full Face Mask - Small
- SFF1002 - Siesta Full Face Mask - Medium
- SFF1003 - Siesta Full Face Mask - Large

## Replacement Supplies Part Numbers

- SFF3001 - Replacement Cushion - Small
- SFF3002 - Replacement Cushion - Medium
- SFF3003 - Replacement Cushion - Large
- SFF2001 - Replacement Headgear
- SFF2002 - Replacement Frame



REACTHEALTH

[www.reacthealth.com](http://www.reacthealth.com)  
[info@reacthealth.com](mailto:info@reacthealth.com)  
Phone: 863-226-6285  
Fax: 863-226-6284