

CPAP XCHANGE
SLEEP WELL TONIGHT

DeVilbiss[®]
HEALTHCARE



Hybrid
Nasal/Oral
CPAP Mask
(HYB500)

Hybrid™ Universal Interface

A Fusion of Comfort and Innovation

The Hybrid is a unique fusion of technologies, designed to be a **universal interface for patients, clinicians, and providers**. This one product serves as a nasal interface, an oral interface, and a combination of the two. The **Hybrid is a dual – airway interface** for use with positive airway pressure devices for the treatment of respiratory insufficiencies and obstructive sleep apnea.

The unique design of the Hybrid promotes patient comfort by **eliminating pressure points** on the forehead and across the bridge of the nose.

The Hybrid complete system **contains all of the components required to properly size the patient**, eliminating the hassle of opening multiple packages to get the right fit.

Performance Features

- Built-in chin flap **gently supports the chin** and keeps the mouth closed; eliminating the need for a chin strap
- Premium headgear provides five points of adjustability to provide **better fit, seal and increased support**
- Quick connect/disconnect on headgear enables **easy removal and application**
- Patients can wear their glasses while wearing the interface

Design Features

- 100% latex free
- One package contains all of the components required to properly size the patient, eliminating the hassle and waste of opening multiple products to properly fit the patient
- **Universal design simplifies choice and fitting protocols** and minimizes time spent with a patient to address issues such as mouth breathing or claustrophobia
- Adjustable nasal pillow design provides a **customized fit for each patient**
- No nasal bridge or forehead contact eliminates nasal bridge sores and air leaks near eyes
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Supplemental oxygen ports
- Patient Instructional DVD and carrying case included with each complete system
- To view video, please visit: <http://www.hybridmask.com>

Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
HYBRID MASK (HYB500)	Combination oral/nasal mask, used with continuous positive airway pressure device. each	A7027	1 every 3 months	\$179.35-\$186.52
HYBRID ORAL CUSHION (HYB511, HYB513, HYB515)	Oral cushion for combination oral/nasal mask, replacement only. each	A7028	2 pairs per month	\$49.54 -\$49.54
HYBRID NASAL PILLOW (HYB521, HYB523, HYB525)	Nasal pillows for combination oral/nasal, replacement only. pair	A7029	2 pairs per month	\$20.24 -\$20.24
HYBRID HEADGEAR (HYB530)	Headgear	A7035	One per 6 months	\$33.75 - \$33.79

Accessories and Replacement Parts

PART NUMBER	PRODUCT DESCRIPTION	PART NUMBER	PRODUCT DESCRIPTION
HYB500	Hybrid Interface Complete System*	HYB530	Hybrid Interface Headgear
HYB515	Hybrid Interface Oral Cushion Large	HYB525	Hybrid Interface Nasal Pillows Large
HYB513	Hybrid Interface Oral Cushion Medium	HYB523	Hybrid Interface Nasal Pillows Medium
HYB511	Hybrid Interface Oral Cushion Small	HYB521	Hybrid Interface Nasal Pillows Small

* Includes small, medium, large oral cushions, small, medium large nasal pillows and headgear.



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com